

<b>IHE Integration Statement</b>		<b>Date</b>
		08 Mar 2008
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
Kowa Optimed, Inc.	VK2	5.0.0.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Eye Care Workflow	Acquisition Modality	PPS Exception Management
		Eye Care Image
Charge Posting	Acquisition Modality	None
<b>Internet address for Kowa Optimed's IHE Information:</b>		
<a href="http://www.kowa-usa.com/frontend/landing_medical.asp">http://www.kowa-usa.com/frontend/landing_medical.asp</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	Not Applicable	
<b>DICOM</b>	<a href="http://www.kowa-usa.com/frontend/landing_medical.asp">http://www.kowa-usa.com/frontend/landing_medical.asp</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>